



Sifam

Medical

Catálogo



Tubos y sondas

A background image showing medical professionals in a clinical setting. A woman in a white lab coat and blue gloves is working with a medical device. A man in a white lab coat and blue gloves is holding a blue clipboard. The image is dimmed to allow text to be overlaid.


PRODUCTOS SOCIOSANITARIOS DE ALTA CALIDAD

En SIFAM Medical ofrecemos soluciones sociosanitarias de primer nivel que mejoran la calidad de vida de los pacientes y facilitan el trabajo de los profesionales de la salud. Descubre nuestra amplia gama de productos diseñados para cumplir con las más altas exigencias del sector.

TUBOS Y SONDAS


Disponemos de una amplia gama de tubos y sondas médicas, diseñados para garantizar la máxima seguridad y precisión en procedimientos de drenaje, ventilación y administración de fluidos.



A close-up, low-angle shot of medical equipment. In the foreground, a clear plastic drip chamber is connected to a clear plastic tube that hangs down. Above it, a syringe with a green plunger is visible, also connected to the tubing. The background is dark and out of focus, showing more of the medical stand.

TUBOS Y SONDAS MÉDICAS

Suction Catheter -Thumb Control (Cap-cone Type)



Ref. GT017-100B

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes
- (4) Connector(color coded)-Size(Fr/Ch): F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.

Suction Catheter (T Type)



Ref. GT017-100C

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade),
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch): F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.


Suction Catheter (Y Type)



Ref. GT0/7-100D

Especificaciones


- (1) Made of non-toxic, non-irritant PVC (Medical grade).
- (2) Color: Transparent translucent (surface fog form).
- (3) Distal end is opened with two lateral eyes.
- (4) Connector (connector) - Size (Fr/Ch): F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



Options Available:

- More connector: Transparent connector.
- More tip: declined tip, big tip and so on.
- With/without X-ray detectable thread.
- Length: 45cm-55cm or any length meet customers' demands.
- Or DEHP free lube.

Extractor Mucus



Ref. GT0/7-200

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free available.
- (2) Suction tube length: 35-40cms.
- (3) Capacity: 25ml.
- (4) Size: Fr8, Fr10, Fr12, Fr14, also available as per customer request
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.

Closed Suction Catheter



Ref. GT023-200

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade) .
- (2) Size(Adult and pediatric): 6Fr, 8Fr, 10Fr, 12Fr, 14Fr, 16Fr.
- (3) Package: Single blister package.
- (4) Sterilized by EO GAS.

Feeding Tube



Ref. GT0JS-100

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade). or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent (surface log form).
- (3) Closed tip with two lateral eyes.
- (4) With/without X-ray detectable thread.
- (5) Connector(color, coded)-Size(Fr/Ch): F4, F5, F6, F8, F10, F12, F14, F16, F18
- (6) Two types connectors available.
- (7) Package:Single in PE or blister package.
- (8) Sterilized by EO GAS.

PVC Stomach Tube



Ref. GT0I9-J00A

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free.
- (2) Color: Transparent or translucent (surface fog form),
- (3) With/without X-ray detectable thread.
- (4) Connector(color coded)-Size(Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (5) Package:Single in PE or blister package.
- (6) Sterilized by EO GAS.

100% Silicone Stomach Tube




Ref. GT019-100B

Especificaciones

- (1) Made of non-toxic, non-irritant 100% silicone.
- (2) Color: Transparent or translucent (surface log form)
- (3) With/without X-ray detectable thread.
- (4) Connector (color coded)-Size (Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.

Duodenal Tube



Ref. GT019-100C

Especificaciones

- (1) Made of non-toxic, non-irritant PVC (Medical grade), or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent (surface fog form).
- (3) With X-ray detectable thread.
- (4) Closed distal end with four lateral eyes.
- (5) Connector (color coded)-Size (Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (6) Package: Single in PE or blister package.
- (7) Sterilized by EO GAS.

Ryle'sTube




Ref. GT019-100D

Especificaciones

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent(surface fog form).
- (3) With X-ray detectable thread.
- (4) With metal ball.
- (5) Connector(color coded)-
Size(Fr/C): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (6) Package: Single in PE or blister package,
- (7) Sterilized by EO GAS.



Rectal Tube




Ref. GT020-100

Especificaciones

- (1) Made of non-toxic, non-irritant PVC.
- (2) Color: Transparent or translucent.
- (3) Closed lip with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch): F20, F22, F24, F26, F28, F30, F32, F34, F36.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.

Nelaton Catheter



GT021-100



Especificaciones

- (1) Made of non-toxic, non-irritant PVC
- (2) Color: Transparent or translucent
- (3) Closed tip with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (6) Nelaton or tiemann, male or female available.
- (7) Package: Single in PE or blister package.
- (8) Sterilized by EO GAS.



**TU SOCIO DE CONFIANZA EN EL
SECTOR SOCIOSANITARIO.**



 www.sifammedical.com
 **900 447 447**
 **673 301 127**