



Sifam

Medical

Catálogo



Urología


A background image showing medical professionals in a clinical setting. A woman in a white lab coat and blue gloves is working with a medical device. A man in a white lab coat and blue gloves is holding a blue clipboard. The image is dimmed to allow text to be overlaid.

PRODUCTOS SOCIOSANITARIOS DE ALTA CALIDAD

En SIFAM Medical ofrecemos soluciones sociosanitarias de primer nivel que mejoran la calidad de vida de los pacientes y facilitan el trabajo de los profesionales de la salud. Descubre nuestra amplia gama de productos diseñados para cumplir con las más altas exigencias del sector.

UROLOGÍA

Ofrecemos una gama de bolsas de orina de alta calidad, diseñadas para garantizar comodidad, seguridad y eficiencia en el manejo de pacientes con necesidades urológicas en entornos clínicos y hospitalarios.



A close-up photograph of a medical professional, likely a urologist, wearing a white lab coat and white gloves. The doctor is performing a physical examination on a patient's lower back and pelvic region. The patient is lying down, and their light-colored pants are pulled down to the knees. The doctor's hands are positioned on the patient's lower back, one hand on the left and one on the right, near the waistline. The background is slightly blurred, showing a clinical setting.

UROLOGÍA



Ref. GT028-100

Economic Urinary Drainage Bag

Capacity: 2000ml

**Without bottom outlet; with 90cm inlet tube,
with non-return valve,**

Sterilization: EO

Packing: 1pc/PE bag, 250pcs/carton

Measurement of carton: 520x320x360mm



Ref. GT028-100

Economic Urinary Drainage Bag

Capacity: 2000ml

with bottom outlet; with 90cm inlet tube;

with non-return valve

with pull-push valve

Sterilization: EO

Form of packing: 1pc/PE bag, 250pcs/carton



Ref. GT028-210

Economic Urinary Drainage Bag

Capacity: 2000ml

**With bottom outlet; with 90cm inlet tube;
with non-return valve, with T valve**

Sterilization: EO

Packing: 1pc/PE bag, 250pcs/carton

Measurement of carton: 520x370x320mm



Ref. GT028-300

Luxury Urinary Drainage Bag

Capacity: 2000ml

**Type: square bottom, latex free, with pull push valve,
with white sample port, with anti-reflux valve, with plastic tie**

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Packing: 1pc/ PE or blister pack, 100pcs/carton,

Measurement of carton: 540x360x320mm



Ref.GT028-310

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: Square bottom; Latex free

with twist valve; with white sample port;

with anti-reflux valve ; with plastic tie

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 540x360x320mm



Ref. GT028-320

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: Square bottom; Latex free

with T valve; with white sample port;

with anti-reflux valve ; without plastic tie

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 540x360x320mm



Ref. GT028-330

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: Square bottom; Latex free

with T valve; with white sample port; with anti-reflux valve

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 540x360x320mm



Ref. GT028-340

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: Square bottom; Latex free

with T valve; without sample port; with anti-reflux valve

with double ropes

Thickness film: 0.11mm per ply, Front film is transparent

Back film is white

Inlet tube: O.D 6.5mm; 90cm length

Sterilization: EO

Form of packing: 1pc/PE pack, 100pcs/carton

Measurement of carton: 540x360x320mm



Ref. GT028-400

Luxury Urinary Drainage Bag

Capacity: 2000ml

**Type: Square bottom; Latex free
with side valve; with anti-reflux valve
with single plastic hanger and tubing clip**

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x320x380mm



Ref. GT028-410

Luxury Urinary Drainage Bag

Capacity: 2000ml

**Type: Square bottom; Latex free
with Single reversal valve; with needle sample port
with anti-reflux valve;
with double plastic hanger and single rope
with blue bed sheet clip**

Inlet tube: O.D 9mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x320x380mm



Ref. GT028-500

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: Square bottom; Latex free

with T valve; with white sample port; with anti-reflux valve

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 540x360x320mm



Ref. GT028-340

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: semicircular bottom, Latex free

with twist valve; with anti-reflux valve; with plastic tie

Front and back films are transparent

Inlet tube: O.D 6.5mm; 90cm length

Sterilization: EO

Form of packing: 1pc/PE pack, 100pcs/carton

Measurement of carton: 440x360x240mm



Ref. GT028-600

Luxury Urinary Drainage Bag

Capacity: 2000ml

Type: pyriform; latex free, with 100% silicone blue outlet and with folding clamp, with needle sample port, with anti-reflux valve, with finger guard clamp

Inlet tube: O.D 9mm; 100cm length

Sterilization: EO

Form of packing: 1pc/PE or blister pack, 20pcs/carton.

Measurement of carton: 380x220x150mm



Ref. GT028-700

Urinary Drainage Bag

Capacity: 2000ml

**Type: Pyriform; Latex free
with T valve; with needle sample port;
with anti-reflux valve
with double plastic hanger and rope;
with white bed sheet clamp**

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT028-800

Urinary Drainage Bag

Capacity: 2000ml

Type: Pyriform; Latex free

with T valve; with needle sample port;

with anti-reflux valve

with blue plastic hanger and rope;

with white bed sheet clamp

Inlet tube: O.D 9mm; 100cm length

Sterilization: EO

Form of packing: 1pc/PE pack, 20pcs/carton

Measurement of carton: 310x260x220mm



Ref. GT028-810

Urinary Drainage Bag

Capacity: 2000ml

Type: Pyriform; Latex free

with T valve;with needle sample port;

with anti-reflux valve

with white double plastic hanger and rope;

with white bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT028-900

Luxury Urinary Drainage Bag

Capacity: 4000ml

with T valve; Latex free;

with anti-reflux drip chamber (one part)

with needle sample port and tubing clip is optional

with anti-reflux valve;

with reinforced double plastic hanger and rope

with blue bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT028-910

Luxury Urinary Drainage Bag

Capacity: 2000ml

with T valve; Latex free;

with anti-reflux drip chamber (one part)

with needle Sample port and tubing clip is optional

with anti-reflux valve;

with reinforced double plastic hanger and rope

with white bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT028-920

Luxury Urinary Drainage Bag

Capacity: 2000ml

with T valve; Latex free;

with anti-reflux drip chamber (Three part)

with needle sample port and tubing clip is optional

with anti-reflux valve;

with reinforced double plastic hanger and rope

with blue bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT029-100

Luxury Urinary Drainage Bag

Capacity: 2000ml

with T valve; Latex free;

with anti-reflux drip chamber (Three part)

with needle sample port and tubing clip is optional

with anti-reflux valve;

with reinforced double plastic hanger and rope

with blue bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 40pcs/carton

Measurement of carton: 430x300x380mm



Ref. GT029-300

Urine Meter Drainage Bag

Capacity: 2600ml +400ml; Latex free

with T valve on urine bag;

with twist valve on urine meter

with needle sample port; with anti-reflux valve

with reinforced double plastic hanger and rope

with white bed sheet clamp

Inlet tube: O.D 10mm; 100cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 10pcs/carton

Measurement of carton: 480x410x250mm



Ref. GT029-310

Urine Meter Drainage Bag

Capacity: 2600ml +500ml, Latex free

with T valve on urine bag ; with twist valve on urine meter

with needle Sample port; with anti-reflux valve

with reinforced double plastic hanger and rope

with white bed sheet clamp

Inlet tube: O.D 10mm; 130cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 10pcs/carton

Measurement of carton: 480x410x250mm



Ref. GT029-320

Urine Meter Drainage Bag

Capacity: 2000ml +200ml; Latex free

with T valve on urine bag; with twist valve on urine meter

with needle sample port; with anti-reflux valve

with reinforced double plastic hanger and rope

with white bed sheet clamp

Inlet tube: O.D 10mm; 130cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 10pcs/carton

Measurement of carton: 480x410x250mm



Ref. GT029-330

Urine Meter Drainage Bag (Double Chambers)

Easy to read molded plastic meters are calibrated up to

400 ml urine meter

with T valve on urine bag ; with twist valve on urine meter

with needle sample port; with anti-reflux valve

with reinforced double plastic hanger and rope

with white bed sheet clamp

Inlet tube: O.D 10mm; 130cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 10pcs/carton

Measurement of carton: 480x410x250mm



Ref. GT029-340

Urine Meter Drainage Bag

Easy to read molded plastic meters are calibrated up to 500 ml urine meter

with 2600ml Luxury Urinary Drainage bag

with anti-reflux drip chamber (one part)

with T valve on urine bag; with needle sample port

with anti-reflux valve;

with four belts & two belts to be choose

Inlet tube: O.D 10mm; 130cm length

Sterilization: EO

Form of packing: 1pc/Blister pack, 10pcs/carton

Measurement of carton: 580x290x235mm



Ref. GT029-400

Leg Bag

Capacity: 600ml

Type: Square

with blue pull-push valve & blue tip cover

with anti reflux valve; with double high elastic belts

Sterilization: EO

Form of packing: 1pc/PE bag, 48pcs/carton

Measurement of carton:

300x260x150mm With 120cm tube



Ref. GT029-410

Leg Bag

Capacity: 500ml & 1000ml

Type: ellipse shape

with green twist valve & green tip cover

with anti reflux valve; with double high elastic belts

Sterilization: EO

Form of packing: 1pc/PE bag, 48pcs/carton

Measurement of carton: 300x260x150mm



Ref. GT029-420

Leg Bag (3 chambers)

Capacity: 500 ml

Type: ellipse shape

with T valve & short outlet tube;

with anti reflux valve

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 550x400x325mm



Ref. GT029-430

Leg Bag

Capacity: 500 ml

Type: ellipse shape
with T valve & short outlet tube;
with double high elastic belts

Sterilization: EO

Form of packing: 1pc/Blister pack, 100pcs/carton

Measurement of carton: 380x260x230mm



Ref. GT029-440

Leg Bag

Capacity: 350/500/750 ml

Type: arc shape
with T valve & lever tap short outlet tube;
with double high elastic belts
Tube length 5~90cm for customer option

Sterilization: EO

Form of packing: 1pc/PE or Blister pack, 200pcs/carton

Measurement of carton: 460x440x300mm

A: SIZE 1.8cm/wide, 48cm/length X 2 Material wo/buttons, Polyester, Nylon

B: SIZE 2.3cm/wide, 50cm/length X 2 Material wo/velcro, Polyester, Nylon



Ref. GT029-450

Leg Bag (3 chambers)

Capacity: 350/500/750 ml

Type: arc shape

with T valve & lever tap short outlet tube;

with double high elastic belts

Tube length 5~90cm for customer option

Sterilization: EO

Form of packing: 1pc/PE or Blister pack, 200pcs/carton

Measurement of carton: 460x440x300mm

A: SIZE 1.8cm/wide, 48cm/length X 2 Material wo/buttons, Polyester, Nylon

B: SIZE 2.3cm/wide, 50cm/length X 2 Material wo/velcro, Polyester, Nylon



Ref. GT029-460

Leg Bag

(2 chambers/3 chambers)

Capacity: 350/500/750 ml

Type: diamond shape

with T valve & lever tap short outlet tube;

with double high elastic belts

Tube length 5~90cm for customer option

Sterilization: EO

Form of packing: 1pc/PE or Blister pack, 200pcs/carton

Measurement of carton: 460x440x300mm

A: SIZE 1.8cm/wide, 48cm/length X 2 Material wo/buttons, Polyester, Nylon

B: SIZE 2.3cm/wide, 50cm/length X 2 Material wo/velcro, Polyester, Nylon



Ref. GT029-500

Disposable Paediatric Urine Collector

Size: 100ml, 200ml

Component: pectic paper, PVC bag, sponge

Sterilization: EO

Form of packing: 1pc/PE pack, 5000pcs/carton

Measurement of carton: 550x540x430mm



Ref. GT029-610

Wound Drainage System

Parts: Container, two to the connector, drainage pipe, connecting pipe, needle, non-return valve, and so on.

Main raw materials: PVC and/or silicone rubber According to drainage pipes and containers of different

materials used can be divided into PP, PS,

SS three types According to the capacity of different

containers can be divided into 200ml, 400ml, 500ml and 800ml.

This product is used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.

200ml 20PCS/43x32x19cm, 400ml 12PCS/43x32x19cm

500ml 24PCS/55x32x48cm, 800ml 12PCS/45x32x48cm



Ref. GT029-620

Wound Drainage System


Parts: Container, two to the connector, drainage pipe, connecting pipe, needle, non-return valve, and so on.

Main raw materials: PVC and/or silicone rubber. According to drainage pipes and containers of different materials used can be divided into PP, PS, SS three types.

According to the capacity of different containers can be divided into 200ml, 400ml, 500ml and 800ml. This product is used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.

200ml 20PCS/43x32x19cm, 400ml 12PCS/43x32x19cm

500ml 24PCS/55x32x48cm, 800ml 12PCS/45x32x48cm



Ref. GT029-800

Disposable Colostomy Bag

Size: 20cm x 14cm, 25cm x 15cm

Adhesive area: 10cm x 9cm, 12cm x 12cm

Center opening: 3cm

Sterilization: EO

Expiry date: 3 years

Form of packing: 20pc/PE bag, 5 bags/box, 20boxes/carton

Measurement of carton: 470x330x240mm(20cm x 14cm)

580x350x240mm(25cm x 15cm)



Ref. GT029-900

Disposable Enema Bag

Capacity: 1500ml

Material: PVC

Sterilization: EO

Expiry date: 3 years

Form of packing: 1pc/PE bag, 200pcs/carton

Measurement of carton: 460x440x330mm



Ref. GT029-900V

Disposable Enema Bag

Capacity: 1500ml

Material: PVC

Sterilization: EO

Form of packing: 1pc/PE bag, 200pcs/carton

Measurement of carton: 460x440x330mm



GT030-200
Simple Type



GT030-300
Foam Type



GT030-400
Glue Pan Type



GT030-500
Two-piece Type



GT030-600
Foam Type



GT030-700
Two-piece Type

NEW STYLE OSTOMY BAG

Our ostomy bag are of excellent quality and advanced usage, making the wound healed in the good condition. Ostomy bag applications: collect excrement and urine after colostomy, ileostomy and urostomy.

Features:

1.The hydrocolloid skin barrier has the performance of anti-maceration, low allergy, and warping and leakage. It has excellent exudation-absorbent and long sticking ability(up to 7days).And it is soft and comfortable.

2.Hydrocolloid skin barrier combines with the highly permeable PU film makes it stick for a long time, no leakage of exudates, convenient to use.

3.The material of ostomy bag possesses high obstructing ability, it is comfortable and covert.

4.The one-piece system design discharge intestinal gas easily and filtering effect it excellent.

Simple type
Foam type
Glue Pan type
Two-piece type

Packing Reference:
10pcs/bag,400pcs/ctn
Size: 350x310x400mm
G.W: 6kg N.W: 4.5kg



**TU SOCIO DE CONFIANZA EN EL
SECTOR SOCIOSANITARIO.**



 www.sifammedical.com
 **900 447 447**
 **673 301 127**